EyeKon Medical, Inc.

2451 Enterprise Road, Clearwater, Florida 33763

ENGLISH HYDROPHILIC ACRYLIC INTRAOCULAR LENSES NATURAL CLEAR/NATURAL YELLOW

Intended Use

The ultra-violet absorbing co-aspheric intraocular lenses manufactured by EyeKon Medical, Inc. are optical implants for the replacement of the human lens in the visual correction of aphakia. The intraocular lens is manufactured from a medical grade co-polymer of Hydrophilic Acrylic, with a polymerisable UV blocker. The lens is designed for implantation in the capsular bag following extracapsular cataract extraction. The optical portion of the lens is designed for folding prior to insertion, allowing the lens to be inserted through an incision size between 3.2 and 3.8 millimeters (mm) when using insertion forceps, or as low as 1.8mm depending on the lens inserter being used. The labeled dioptic power of the lens is in aqueous (balance solutions).

How the Device is Supplied

Each intraocular lens is supplied in a sterile tamperproof borosilicate glass vial and sealed within a sterile pouch.


Indications

The intraocular lens is designed to be used for visual correction of aphakia for primary implantations following extracapsular cataract extraction. Implantation of the intraocular lens should not be performed in patients under 18 years of age.


Spectral Transmission Curve

The cutoff wavelength of 10% transmission for UV light is 380 nm for +20.0 diopter lens.

Spectral Transmission Curve SC25-FOLD UV

Spectral Transmission Natural Yellow Hydrophilic Acrylic Lens

Instructions For Use

- Prior to implanting, examine the lens package for type, power and proper configuration.
- In a sterile environment, peel open the lens pouch to attain the sterile vial.
- To remove screw cap, twist counterclockwise.
- To avoid dehydration, leave lens immersed in vial until ready to fold and implant.
- Remove the lens from vial.
- Place the lens in a petri dish with BSS solution. The anterior side of the lens should be facing up. Grasp the lens at the 3 and 9 o'clock position with the folder. Fold the lens by squeezing the folding forceps while keeping the lens against the bottom of the petri dish.
- Thoroughly rinse the lens with BSS solution.
- During insertion, the haptic tip should point to the surgeon's left and should enter the incision first. The inferior haptic should first be placed beneath the capsular rim. The forceps should be rotated clockwise and opened slowly. The superior haptic will rotate outside of the eye. Once the lens is unfolded, introduce the superior haptic into the capsular bag in a fashion similar to a conventional lens.
- If using a lens injector, follow manufacturer's instructions.

Contraindications

Patients with certain conditions are not suitable candidates to receive a hydrophilic acrylic intraocular lens because the lens may exacerbate an existing condition, or may cause an unreasonable risk to the patient's eyesight. They include the following:

- Prior intraocular surgery in the operative eye.
- Iris atrophy
- Proliferative diabetic retinopathy
- Chronic, medically uncontrolled glaucoma.
- Chronic use of steroids, immunosuppressive and/or antineoplastic agents.
- Chronic uveitis, iritis, iridocyclitis, or rubeosis of the lids.

Cautions

As with any surgical procedure, there is a potential risk involved. Potential complications accompanying intraocular lens implantations may include, but are not limited to the following:

Wound Leak	Corneal Edema	Retinal Detachment
Pupillary Block	Cystoid Macular Edema	Vitritis
Decentration of the Implant	Corneal Decompensation	
Endophthalmitis	Dislocation and Subluxation of the Implant	

Some of the listed complications may require secondary surgical intervention. The complications could be observed following implantation of any intraocular lenses. A temporary (30-60 minutes) opacification of the lens may occur if the temperature of the lens is below 25_ C prior to implantation. This opacification occurs as the lens equilibrates from a low temperature to body temperature. This may be minimized by allowing the lens to reach 25°C before implanting.

Warnings

The safety of intraocular lens implantation has not been substantiated in patients with pre existing ocular conditions (chronic drug miosis, glaucoma, amblyopia, diabetic retinopathy, previous corneal transplant, history of retinal detachment, etc.) Physician considering lens implantation in such patients should explore the use of alternative methods of aphakic correction. The long term effects of intraocular lens implantation have not been determined. Therefore, physicians should continue to monitor implant patients postoperatively on a regular basis. The safety and effectiveness of this lens, if placed in the anterior chamber, has not been established. Implantation of posterior chamber lenses in the anterior chamber has been shown in some cases to be unsafe. Improper handling or folding techniques may cause damage to the haptic or optic portions of the lens. Physicians should not attempt to implant lenses that have optical tears or broken haptics. The aqueous in the natural eye ranges from 7.0 to 7.5 in pH. U.S. Pharmacopoeia has required that all pH balance range from 6.8 to 7.6. The high pH balance that can be found in certain viscoelastic, balance salt solutions and pharmaceutical drugs can affect the pH balance of the hydrophilic acrylic lens. EyeKon Medical, Inc. would recommend that physicians do not deviate from 6.8 to 7.6 in all BSS, viscoelastics and all ophthalmic eye drops including steroids. Changes in the clarity of the hydrophilic acrylic lens could occur should the pH balance differ (i.e. lens could become yellow or brown).

Precautions

- Do not resterilize the lens. The sterility of the lens cannot be assured unless sterilized by EyeKon Medical. Lenses requiring resterilization should be returned to the manufacturer.
- Do not use the device if the sterile packaging has been damaged as sterility may have been compromised.
- Do not soak lenses in solutions other than balanced saline or equivalent as fluids may be absorbed by the lens.
- Once the packaging is no longer intact the intraocular lens must be used without delay to avoid dehydration of the lens.
- Because of its hydrophilic nature, the lens may absorb substances with which it comes into contact (disinfectant, medicines, blood cells). This may cause a "Toxic Lens Syndrome". Rinse the lens carefully prior to implantation into the eye.
- It is imperative that the hydrophilic acrylic lens is implanted in the capsular bag.
- A high level of surgical skill is required for intraocular lens implantation. A surgeon should have observed and/or assisted in numerous surgical implantations. Do not use after the expiration date on the package. The surgeon must be aware of the risk of opacification of the intraocular lens, which may necessitate lens removal. Any case involving lens removal must be reported to EyeKon Medical, Inc.

Adverse Reactions

The following adverse reactions have been reported following cataract extraction and implantation of an intraocular lens:
Phophopyn, Intraocular Infection, Corneal Decompensation, Lens Removal, Pupillary Block and Lens Dislocation and/or Replacement.

Calculation of Lens Power

The power of the implanted intraocular lens for each patient can be estimated from the earlier refractive error or calculated from the corneal radius, the depth of the anterior chamber and the axial length of the eye.

Return of Goods and Resterilization Policy

Unopened lenses will be exchanged for a comparable dollar value provided they have not exceeded their expiration date. The returned lens must be in its original, unopened, and undamaged packaging to be eligible for replacement or credit. A reprocessing fee may be assessed for the unopened lenses. Prior to returning the

lenses the customer must obtain a Return Authorization number, and provide serial number, power, and reason for return. It is EyeKon's policy not to issue credit or cash refunds for returned lenses.

DO NOT ATTEMPT TO RESTERILIZE THE LENSES.

Reporting

Adverse reactions and/or potentially sight threatening complications that may reasonably be regarded as lens related and that were not previously expected in nature, severity or degree of incidence should be reported to EyeKon Medical, Inc., 2451 Enterprise Road, Clearwater, Florida 33763 at 727-793-0170 or 1-800-633-9248 domestically.

This information is being requested from all implant surgeons. Physicians are encouraged to report these events in order to aid in identifying emerging or potential problems with intraocular lenses. These problems may be related to a specific lot of lenses or may be indicative of long-term effects associated with these lenses or with IOLs in general.

Disclaimer of Liability

EyeKon will not be liable for any injury suffered to a patient as a result of: Any implantation method or technique used by a physician to implant the lens. Any prescription selection and use of the lens for any individual patient or patient's condition. EyeKon makes no expressed or implied warranties in connection with the sale of this intraocular lens and specifically disclaims any warranty liability of the marketability or fitness for use. CAUTION: Federal (USA) law restricts this device to sale by or on the order of physician.

Bibliography

• Allarakhia, L. R.L. Knoll and R.L. Lindstrom. "Soft Intraocular Lenses" Cataract Refract Surge. 113(11): 607-620, 1987. •Utrata, P.J. "Implantation of STAAR Three-Piece Silicone Elastimide Lens." Small-Incision Cataract Surgery. Gills and Sanders ed. 1990 SLACK Inc. Thoroface, NJ 89-101. •Popham, J.K., D.J. Apple, D.A. Newman, R.A. Isenberg and J. Deacon "Advantages and Limitations of Soft Intraocular Lenses: A Scientific Perspective." Soft Implant Lenses in Cataract Surgery. Mazzocco, Rajachik and Epstein ed.1986 Slack In. Thorofare, NJ. 11-30. •Jaffe, N.S. "Current concepts in posterior chamber lens technology." American Intraocular Implant Society Journal, 1985, 11, 456-460. •Binkhorst, C.D., Manschot, W.A. "The mechanism of capsular fixation of intraocular lenses in the human eye." American Intraocular Implant Society Journal, 1977, 3, 100-115. •Apple, D.J. et al. "Loop fixation of posterior chamber intraocular lenses." Cataract, 1984, 7-10.

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Bibliographie

• Allarakhia, L. R.L. Knoll and R.L. Undstrom. "Soft Intraocular Lenses" Cataract Refract Surge. 113(11): 607-620, 1987. •Utrata, P.J. "Implantation of STAAR Three-Piece Silocone Elastimide Lens." Small-Incision Cataract Surgery. Gills and Sanders ed. 1990 SLACK Inc. Thoroface, NJ 89-101. •Popham, J.K., O.J. Apple, D.A. Newman, R.A. Isenberg and J. Deacon "Advantages and Limitations of Soft Intraocular Lenses: A Scientific Perspective." Soft Implant Lenses in Cataract Surgery. Mazzocco, Rajachik and Epstein ed.1986 Slack In. Thorofare, NJ. 11-30. •Jaffe, N.S. "Current concepts in posterior chamber lens technology." American Intraocular Implant Society Journal, 1985, 11, 456-460. •Binkhorst, C.D., Manschot, W.A. "The mechanism of capsular fixation of intraocular lenses in the human eye." American Intraocular Implant Society Journal, 1977, 3, 100-115. •Apple, D.J. et al. "Loop fixation of posterior chamber intraocular lenses." Cataract, 1984, 7-10Calcul de la puissance
•Binkhorst, R.D., Intraocular Lens Power Calculation Manual. New York: Richard D. Binkhorst, 1978. •Donzis, P.B., et al. "An intraocular lens formula for short, normal and long eyes." CLAO Journal, 1985, 11 (2) 95-98. •Gills, J.P. "Minimizing postoperative refractive error." Contact Lens, 1980, 6 (1), 56-59. •Hoffer, K.J. "Preoperative evaluation of the cataractous patient." Survey of Ophthalmology, 1984, 29 (1) 55-69. •Holladay, J.T., et al. "Improving the predictability of intraocular lens power calculation." Arch. Ophthalmology, 1986, 104, 538-541. •Retzlaf, J. et al. "A Manual of Implant Power Calculation SRK Formula, 1981."Richards, S.C. et al. "Clinical evaluation of six intraocular lens calculation formulas." 1985, 11, 153-158.

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Literatur

• Allarakhia, L. R.L. Knoll and R.L. Lindstrom. "Soft Intraocular Lenses" Cataract Refract Surge. 113(11): 607-620, 1987.
•Ultrata, P.J. "Implantation of STAAAR Three-Piece Silicone Elastimide Lns." Small-Incision Cataract Surgery. Gills and Sanders ed. 1990 SLACK Inc. Thoroface, NJ 89- 101.
•Popham, J.K., D.J. Apple, D.A. Newman, R.A. Isenberg and J. Deacon "Advantages and Limitations of Soft Intraocular Lenses: A Scientific Perspective." Soft Implant Lenses in Cataract Surgery. Mazzocco, Rajacich and Epstein ed.1986 Slack In. Thorofare, NJ. 11-30.
•Jaffe, N.S. "Current concepts in posterior chamber lens technology." American Intraocular Implant Society Journal, 1985, 11, 456-460.
•Binkhorst, C.D., Manschot, W.A. "The mechanism of capsular fixation of intraocular lenses in the human eye."American Intraocular Implant Society Journal, 1977, 3, 100-115.
•Apple, D.J. et al. "Loop fixation of posterior chamger intraocular lenses." Cataract, 1984, 7-10.

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• Binkhorst, R.D., Intraocular Lens Power Calculation Manual. New York: Richard D. Binkhorst, 1978.
•Donzis, P.B., et al. "An intraocular lens formula for short, normal and long eyes." CLAO Journal, 1985, 11 (2) 95-98.
•Gills, J.P. "Minimizing postoperative refractive error." Contacta Lens, 1980, 6 (1), 56- 59.
•Hoffer, K.J. "Preoperative evaluation of the cataractous patient." Survey of Ophthalmology, 1984, 29 (1) 55-69.
•Holladay, J.T., et al. "Improving the predictability of intraocular lens power calculation." Arch. Ophthalmology, 1986, 104, 538-541.
•Retzlaf, J. et al. "A Manual of Implant Power Calculation SRK Formula, 1981.
• Richards, S.C. et al. "Clinical evaluation of six intraocular lens calculation formulas." 1985, 11, 153-158

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• Allarakhia, L. R.L. Knoll and R.L. Lindstrom. "Soft Intraocular Lenses" Cataract Refract Surge. 113(11): 607-620, 1987.
•Ultrata, P.J. "Implantation of STAAAR Three-Piece Silicone Elastimide Lens." Small-Incision Cataract Surgery. Gills and Sanders ed. 1990 SLACK Inc. Thoroface, NJ 89- 101.
•Popham, J.K., D.J. Apple, D.A. Newman, R.A. Isenberg and J. Deacon "Advantages and Limitations of Soft Intraocular Lenses: A Scientific Perspective." Soft Implant Lenses in Cataract Surgery. Mazzocco, Rajacich and Epstein ed.1986 Slack In. Thorofare, NJ. 11-30.
•Jaffe, N.S. "Current concepts in posterior chamber lens technology." American Intraocular Implant Society Journal, 1985, 11, 456-460.
•Binkhorst, C.D., Manschot, W.A. "The mechanism of capsular fixation of intraocular lenses in the human eye."American Intraocular Implant Society Journal, 1977, 3, 100-115.
•Apple, D.J. et al. "Loop fixation of posterior chamger intraocular lenses." Cataract, 1984, 7-10.

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• Binkhorst, R.D., Intraocular Lens Power Calculation Manual. New York: Richard D. Binkhorst, 1978.
•Donzis, P.B., et al. "An intraocular lens formula for short, normal and long eyes." CLAO Journal, 1985, 11 (2) 95-98.
•Gills, J.P. "Minimizing postoperative refractive error." Contact Lens, 1980, 6 (1), 56- 59.
•Hoffer, K.J. "Preoperative evaluation of the cataractous patient." Survey of Ophthalmology, 1984, 29 (1) 55-69.
•Holladay, J.T., et al. "Improving the predictability of intraocular lens power calculation." Arch. Ophthalmology, 1986, 104, 538-541.
•Retzlaf, J. et al. "A Manual of Implant Power Calculation SRK Formula, 1981.
• Richards, S.C. et al. "Clinical evaluation of six intraocular lens calculation formulas." 1985, 11, 153-158.

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

Bibliografia

• Allarakhia, L. R.L. Knoll and R.L. Lindstrom. "Soft Intraocular Lenses" Cataract Refract Surge. 113(11): 607-620, 1987.
•Ultrata, P.J. "Implantation of STAAAR Three-Piece Silicone Elastimide Lns." Small-Incision Cataract Surgery. Gills and Sanders ed. 1990 SLACK Inc. Thoroface, NJ 89-101.
•Popham, J.K., D.J. Apple, D.A. Newman, R.A. Isenberg and J. Deacon "Advantages and Limitations of Soft Intraocular Lenses: A Scientific Perspective." Soft Implant Lenses in Cataract Surgery. Mazzocco, Rajacich and Epstein ed.1986 Slack In. Thorofare, NJ. 11-30.
•Jaffe, N.S. "Current concepts in posterior chamber lens technology." American Intraocular Implant Society Journal, 1985, 11, 456-460.
•Binkhorst, C.D., Manschot, W.A. "The mechanism of capsular fixation of intraocular lenses in the human eye." American Intraocular Implant Society Journal, 1977, 3, 100-115.
•Apple, D.J. et al. "Loop fixation of posterior chamger intraocular lenses." Cataract, 1984, 7-10.

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• Binkhorst, R.D., Intraocular Lens Power Calculation Manual. New York: Richard D. Binkhorst, 1978.
•Donzis, P.B., et al. "An intraocular lens formula for short, normal and long eyes." CLAO Journal, 1985, 11 (2) 95-98.
•Gills, J.P. "Minimizing postoperative refractive error." Contact Lens, 1980, 6 (1), 56-59.
•Hoffer, K.J. "Preoperative evaluation of the cataractous patient." Survey of Ophthalmology, 1984, 29 (1) 55-69.
•Holladay, J.T., et al. "Improving the predictability of intraocular lens power calculation." Arch. Ophthalmology, 1986, 104, 538-541.
•Retzlaf, J. et al. "A Manual of Implant Power Calculation SRK Formula, 1981.
• Richards, S.C. et al. "Clinical evaluation of six intraocular lens calculation formulas." 1985, 11, 153-158.

Symbol	Description
	Attention: See instructions for use
	Do not reuse
	Batch Code
	Sterilization by heat
	Use until

Manufactured By 

EyeKon Medical, Inc.
2451 Enterprise Road
Clearwater, FL 33763
Telephone (727) 793-0170
Fax (727) 799-2212
www.eyekonmedical.com

EyeKon Medical, Ltd.
1 Friary Temple Quay Bristol,
BS1 6EA United Kingdom