## **EyeKon Medical, Inc.**

## Viscoelastic hydroxypropyl methylcellulose (HPMC)

### **ENGLISH**

#### Description

HPMC Viscoelastic Ophthalmic Viscosurgical Device is a clear, sterile, isotonic, non pyrogenic, and non-inflammatory viscoelastic solution based on 2% hydroxypropyl methylcellulose (HPMC) with a nominal viscosity of 3,000 cPs. It is non-organic; water soluble; enzyme resistant, which helps provide viscosity stability during storage; and pH stable. HPMC viscoelastics have been demonstrated to have a high degree of biocompatibility. Its dispersive characteristics are especially suited for protecting the endothellium, while it effectively maintains space and gently manipulates delicate intraocular tissue.

HPMC Concentration:	2.0%
Viscosity:	HPMC- 3, 000 ± 1,500 cPs HPMC PLUS- 16, 000 ± 5, 000 cPs
Molecular Weight:	2 million Daltons
pH:	6.8-7.6
Osmolarity:	200 – 400 mOsM
Particulates:	≤10,000 particles >10µ ≤1,000 particles >25µ
Bacterial Endotoxins:	Not more than 0.5 EU/ml
Sterility: Sterile using method described in USP<71>	

#### Intended use

HPMC Viscoelastic is indicated for use as an ophthalmic surgical aid in the anterior segment in patients undergoing cataract extraction and intraocular lens (IOL) implantation. HPMC Viscoelastic is intended to create and maintain space to manipulate intraocular tissue in the anterior segment during cataract extraction and intraocular lens implantation.

Ophthalmic viscosurgical devices serves to maintain a deep anterior chamber during anterior segment surgery allowing reduced trauma to the corneal endothelium and surrounding ocular tissues. The viscoelastic properties help to push back the vitreous face and prevent formation of a flat chamber postoperatively.

#### Contraindications

At present there are no known Contraindications to the use of Viscoelastic when used as recommended; care should be used in patients with hypersensitivity to any components in this material.

#### Precautions

- Precautions normally associated with anterior segment surgical procedures should be observed.
- It is recommended that Viscoelastic be removed from the eye as is possible by thoroughly irrigating and aspirating with a balanced salt solution at the close of surgery.

#### Warnings

- Do not reuse cannula.
- Use only if material is clear.
- Avoid trapping air bubbles
- Use only if container is undamaged. Contents are only sterile if the pouch is undamaged.
- Do not resterilize the syringe.
- Do not use after expiration date shown on syringe as sterility cannot be assured after the expiration date.

#### Adverse reactions

- HPMC Viscoelastic is tolerated after injection into human eyes during ophthalmic surgical procedures. As with most viscoelastics, a transit rise in intraocular pressure has been reported in some cases.
- Postoperative inflammatory reactions such as hypopyon and iritis have been reported with the use of ophthalmic viscoelastics, as well as incidents of corneal edema and corneal decompensation. Their relationship to the use of HPMC Viscoelastic has not been established.

#### Applications

Cataract Surgery – IOL Implantation

Use sterile opening technique when removing contents of the peel pouches. In a sterile environment, the pouch is peeled open to attain the sterile syringe. A cannula is used to slowly and carefully inject a sufficient amount of Viscoelastic into the anterior chamber. The injection may be performed before or after removal of the crystalline lens. Additional viscoelastic will be injected as needed during anterior segment surgery to fully maintain anterior chamber depth, to manipulate intraocular tissue, or to replace any viscoelastic lost during the surgical procedure. Viscoelastic may also be used to coat surgical instruments and the intraocular lens prior to implantation.

#### How supplied

Viscoelastic is filled into a sterile, single use, pre lubricated glass syringe with Luer lock. A gray butyl stopper is inserted into the syringe, followed by a polycarbonate plunger rod. The assembled syringe is sealed into a peel pouch, steam sterilized and packaged into a box. The complete assembled box is shrink wrapped for protection. Store at approximately 25 °Celsius. Avoid exposure to direct light. Available in 1.0ml-2.5ml, accompanied with a sterile cannula.

#### Directions For Use

##### FOR INTRAOCULAR USE. FOR SINGLE USE ONLY.

- Peel open pouch from both syringe and cannula under aseptic conditions.
- Remove rubber tip cap from syringe tip.
- Thread the cannula onto the syringe by holding the Luer Lock adapter between your thumb and forefinger. Twist the cannula in a continuous motion using the cartridge as a wrench with the other hand. Twist until the cannula hub has traveled the full length of the threads and is firmly seated.
- Visually inspect that the cannula hub has traveled the full length of the Luer Lock adapter.
- Remove the plastic cartridge from the cannula in a straight motion, being sure not to twist or unscrew the cannula while removing the cartridge.
- Purge the remaining air from the system by holding the syringe barrel with one hand and gently depressing the plunger rod with the other until Viscoelastic appears at the cannula tip.

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Made in the U.S.A.
EyeKon Medical Ltd.