EyeKon Medical, Inc.

2451 Enterprise Road, Clearwater, Florida 33763

ENGLISH

Capsular Tension Rings

Introduction

Capsular Tension Rings are manufactured from clinical quality polymethyl methacrylate (PMMA). They are used to stabilize the capsular bag in various situations involving complications during IOL surgery, including:

- damaged or missing zonules
- complications during Phaco surgery
- repositioning of decentered IOLs
- prevention of excessive capsular shrinking (high myopia)
- prevention of IOL subluxation

To allow the surgeon flexibility in meeting individual patient requirements, capsular tension rings are available in 10.0mm and 12.0mm compressed diameters.

Indications

The capsular tension ring effectively stabilizes the capsular bag in cases of weak or partially absent zonules, reducing the rate of serious complications such as vitreous loss, dislocation of the nucleus, and inability to implant a PC IOL.

Cautions

The general risks involved with the use of a capsular tension ring are those common to all surgical procedures involving cataract extraction and intraocular lens implantation. Due to the indicated use for the capsular tension ring, the patients are already at a higher risk for surgical complications. The risk to these "high-risk" patients is actually reduced by the use of the capsular tension ring.

Warnings

Physicians should not attempt to implant capsular tension rings that are broken or have broken eyelets.

Precautions

Do not resterilize the capsular tension rings. The sterility of the capsular tension rings cannot be assured unless sterilized by EyeKon Medical. Capsular tension rings requiring resterilization should be returned to the manufacturer.

Do not soak the capsular tension rings with any solution other than a sterile balanced salt solution or a sterile normal saline. Other solutions may damage the capsular tension rings.

Do not use the capsular tension rings after the expiration date shown on the outside of the package as sterility cannot be guaranteed.

Handle the capsular tension rings carefully. Rough or excessive handling may damage the capsular tension rings.

A high level of surgical skill is required for capsular tension ring implantation. A surgeon should have observed and/or assisted in numerous surgical implantations and successfully completed one or more courses on capsular tension rings prior to attempting to capsular tension rings.

Adverse Reactions

The following adverse reactions have been reported following cataract extraction and implantation of an intraocular lens: Hypopyon, Intraocular Infection, Corneal decompensation, Lens Removal, Pupillary Block and Lens Replacement. Similar adverse reactions may be associated with intraocular lens implantation with capsular tension rings.

How Supplied

Capsular tension rings are ethylene oxide sterilized in a plastic (K-Resin KR-01) case contained within a sealed pouch. The contents of the pouch are STERILE unless the package is damaged or opened.

Directions for Use

Peel open the pouch and remove the CTR. The capsular tension ring should be inspected for debris or damage. Care must be taken to remove the capsular tension ring without bending it into an unnatural state. Breakage should be eliminated by lifting the capsular tension ring straight upward out of the case. You may also release the capsular tension ring by inverting the case and allowing the capsular tension ring to fall gently onto a sterile surface.

Rinse the capsular tension ring with balanced salt solution prior to implantation into the eye.

Physicians should not attempt to implant capsular tension rings that are broken or have broken eyelets.

The packaging box contains peelable labels which display the serial number, model, trade name, diameter and expiration date of the capsular tension ring. These labels are for convenience in maintaining and reporting records of implanted capsular tension rings during the clinical investigation. They are designed to be affixed to the patient's hospital chart, the physician's chart and the implant notification card. One of these labels should be affixed to the patient identification card contained in this package and given to the patient as a permanent record of the implant.

In the occurrence a CTR needs to be destroyed, proper disposal is as follows:

1. Defacing identifying labels (case packaging) by permanent marking.
2. Cutting CTR in half before disposal.

Surgical Technique

A variety of surgical techniques may be employed during the implantation of a CTR. Therefore, the surgeon is best advised to use the method which his training and judgment dictate to be best for the patient. Only tools and equipment intended for use in cataract surgery should be used to implant the CTR. Below are standard techniques that may be followed:

1. Ensure that the capsular bag is fully inflated with an ophthalmic viscoelastic device (OVD). Proceed as directed in the respective viscoelastic IFU.
2. When using an injection device, use only a device which is specifically designed for use with Capsular Tension Rings. Follow the injector manufacturer's instructions regarding preparation and use.
3. During the implantation of the CTR it is slowly "dialed" into the capsular bag through the corneal incision and capsulorhexis either with a forceps or an injector.
4. The direction of the CTR's insertion ("dialing clockwise or anticlockwise") depends on the location of the zonular damage as the tip of ring during the insertion needs to support the region of the zonular weakness/dehiscence rather than further stretching it.

Reporting

Adverse reactions and / or potentially sight threatening complications that may reasonably be regarded as capsular tension ring related and that were not previously expected in nature, severity or degree of incidence should be reported to EyeKon Medical, Inc., 2451 Enterprise Road, Clearwater, FL. 33763. This information is being requested from all implant surgeons.

• Physicians are encouraged to report these events in order to aid in identifying emerging or potential problems with capsular tension rings. These problems may be related to a specific lot of capsular tension rings or may be indicative of long –term effects associated with these capsular tension rings or with capsular tension rings in general.

• If the patient has an EyeKon Medical capsular tension ring and you wish to report, please call: EyeKon Medical, Inc. (800) 633-9248 Local: (727) 793-0170 Fax: (727) 799-2212.

Return of Goods and Resterilization Policy

Opened or unopened capsular tension rings will be exchanged for a comparable dollar value, provided they have not exceeded their expiration date. A reprocessing charge may be assessed for the capsular tension rings that have exceeded their expiration date. This applies to opened or unopened capsular tension rings. It is EyeKon Medical's policy not to issue credit or cash refunds for returned capsular tension rings.

Disclaimer of Liability

EyeKon Medical, Inc. will not be liable for any injury suffered to a patient as a result of:

- Any implantation method or technique used by a physician to implant the capsular tension ring.
- Any selection and use of the capsular tension ring for any individual patient or patient's condition.

EyeKon Medical, Inc. makes no express or implied warranties in connection with the sale of this capsular tension ring and specifically disclaims any warranty liability of the marketability or fitness for use.

Bibliography

Cionni RJ, Osher RH. Management of profound zonular dialysis or weakness with a new endocapsular ring designed for scleral fixation. *Journal of Cataract & Refractive Surgery* 1998; 24(10):1299-1306. Dietlein TS, Jacobi PC, Konen W, Krieglstein GK. Complications of endocapsular tension ring implantation in a child with Marfan's syndrome. *Journal of Cataract & Refractive Surgery* 2000; 26(6):837-940. New and Emerging Techniques – Surgical Intraocular Tension Ring Faschingr OW, Eckhardt M. Complete capsulorhexis opening occlusion despite capsular tension ring implantation. [see comments]. *Journal of Cataract & Refractive Surgery* 1999; 25(7):1013-1015. Gimbel HV, Sun R, Heston JP. Management of zonular dialysis in phacoemulsification and IOL implantation using the capsular tension ring. *Ophthalmic Surgery & Lasers* 1997; 28(4):273-281. Groessl SA, Anderson CJ. Capsular tension ring in a patient with Weill-Marchesani syndrome. *Journal of Cataract & Refractive Surgery* 1998; 24(8):1164-1165. Kugelberg U, Zetterstrom C, Lundgren B, Norrby S, Syren-Nordqvist S. After-ctaract and ocular growth in newborn rabbit eyes implanted with a capsule tension ring. *Journal of Cataract & Refractive Surgery* 1997; 23 Suppl 1:635-640. Lam DS, Young AL, Leung AT, Rao SK, Fan DS, Ng JS. Scleral fixation of a capsular tension ring for severe ectopia lentis. *Journal of Cataract & Refractive Surgery* 2000; 26(4):609-612. Lang Y, Fineberg E, Garzozzi HJ. Vitrectomy to remove a posteriorly dislocated endocapsular tension ring. *Journal of Cataract & Refractive Surgery* 2001; 27(3):474-476.

Lee DH, Lee HY, Lee KH, Chung KH, Joo CK. Effect of a capsular tension ring on the shape of the capsular bag and opening and the intraocular lens. *Journal of Cataract & Refractive Surgery* 2001; 27(3):452-456. Legler U, Witschel B, Lim S, et al. The Capsular Ring: A New Device for Complicated Cataract Surgery. Symposium on Cataract, IOL and Refractive Surgery, Seattle, Washington, USA, 1993.

Menapace R, Findl O, Georgopoulos M, Rainer G, Vass C, Schmetterer K. The capsular tension ring: designs, applications, and techniques. *Journal of Cataract & Refractive Surgery* 2000; 26(6):898-912.

Morcher GmbH, Stuttgart, Germany. http://www.morcher.com . Accessed August 2001. Prevention of Blindness and Deafness Group WHO. A Global Initiative for the Elimination of Avoidable Blindness. http://www.who.int/pbd/Vision2020/V2020slides/index.htm . Accessed August 2001.

Strenn K, Menapace R, Vass C. Capsular bag shrinkage after implantation of an open-loop silicone lens and a poly(methyl methacrylate) capsule tension ring. *Journal of Cataract & Refractive Surgery* 1997; 23(10):1543-1547.

Sun R, Gimbel HV. In vitro evaluation of the efficacy of the capsular tension ring for managing zonular dialysis in cataract surgery. *Ophthalmic Surgery & Lasers* 1998; 29(6):502-505.

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



Bibliografia

Cionni RJ, Osher RH. Management of profound zonular dialysis or weakness with a new endocapsular ring designed for scleral fixation. Journal of Cataract and Refractive Surgery 24(10):1299-1306.
 Dieltien TS, Jacobi PC, Konen W, Krieglsein GK. Complications of endocapsular tension ring implantation in a child with Marfan's syndrome. Journal of Cataract and Refractive Surgery 2000;26 (6):937-940.New and Emerging Techniques –Surgical Intracapsular Tension Ring
 Feschinger CW Eckhardt M. Complete capsulorhexis opening occlusion despite capsular tension ring implantation. [see comments]. Journal of Cataract and Refractive Surgery 1999; 25(7):1013-1015.
 Gimbel HV, Sun R, Heston JP. Management of zonular dialysis in phacoemulsification and IOL implantation using the capsular tension ring. Ophthalmic Surgery and Lasers 1999; 28(4):273-281
 Grossi SA, Anderson CJ. Capsular tension ring in a patient with Weill, Marchesani syndrome. Journal of Cataract & Refractive Surgery 1998; 24(8):1164-1165,
 Kugelberg U, Zetterstrom C, Lundgren B, Norrby S,Syren-Nordqvist S. After Cataract and ocular growth in new-born rabbits eyes implanted with a capsula tension ring. Journal of Cataract and Refractive Surgery 1997;23 Suppl:635-640.
 Lam DS, Young AL, Leung AT, Rao SK Fan OS, Ng JS. Scleral fixation of a capsular tension ring for severe ectopia lentis. Journal of Cataract and Refractive Surgery 2000; 26(4) :609-612.
 Lang Y, Fineberg E, Garzozzi HJ. Vitrectomy to remove a posteriorly dislocated endocapsular tension ring. Journal of Cataract and Refractive Surgery 2001; 27(3):474-476.
 Lee DK, Lee HY, Lee KH, Chung KH, Joo CK. Effect of a capsular tension ring on the shape of the capsular bag and opening and the intraocular lens. Journal of Cataract and Refractive Surgery 2001; 27(3):452-456.
 Legler U, Witschel B, Lim S, et al. The Capsular Ring: A New Device for Complicated Cataract Surgery. Symposium on Cataract, JOL and Refractive Surgery. Seattle, Washington, USA. 1993
 Menapace R, Findl O, Georgopoulos M, Rainer G,Vass C, Schmetterer K The capsular tension ring: designs, applications, and techniques. Journal of Cataract and Refractive Surgery 2000; 26(6):898-912
 Marcher GmbH. Stuttgart, Germany. http://www.morcher.com Accessed August 2001
 Prevention of Blindness and Deafness Group WHO. A Global Initiative for the Elimination of Avoidable Blindness. http://www.who.int/bpd/Vision2020/V2020slides/index.htm. Accessed August 2001.
 Strenn K, Menapace R, Vass C. Capsular bag shrinkage after implantation implantation of an open-top silicone lens and a polymethyl methacrylate capsula tension ring. Journal of Cataract and Refractive Surgery 1997; 23(10): 1543-1547.
 Sun R, Gimbel MV. In vitro evaluation of the efficacy of the capsular tension ring for managing zonular dialysis in cataract surgery. Ophthalmic Surgery and Lasers 1998; 29(6):502-505

Federal (USA) law restricts this device to sale by or on the order of a physician.

Symbol	Description
	Attention: See instructions for use
	Do not reuse
	Batch Code
	Sterilization by EO
	Use until

Manufactured By:

EyeKon Medical, Inc.
 2451 Enterprise Road
 Clearwater, FL 33763
 Telephone (727) 793-0170
 Fax (727) 799-2212

www.eyekonmedical.com

email: info@eyekonmedical.com

Made in the U.S.A.

EyeKon Medical Ltd.
 1Friary Temple Quay
 Bristol, BS1 6EA
 United Kingdom