**EyeKon Medical, Inc.**  
**Viscoelastic Sodium Hyaluronate (SH)**

**ENGLISH**

**Description**

Viscoelastic SH Ophthalmic Viscosurgical Device is a clear, sterile, isotonic, non pyrogenic, and non-inflammatory viscoelastic solution based on sodium hyaluronate (hyaluronic acid). Sodium Hyaluronate is a naturally occurring, high molecular mass, and functions as a tissue lubricant. Viscoelastic SH is suitable for aqueous and vitreous humor and coats the corneal endothelium in ophthalmic surgery. Elasticity of the solution assists in absorbing mechanical stress and providing a protective buffer for tissues. This viscoelasticity enables maintenance of a deep chamber during surgical manipulation since the solution does not flow out of the open anterior chamber. Sodium hyaluronate viscoelastics have been demonstrated to have a high degree of biocompatibility.

HA Concentration: 1.0%, 1.4%, 1.6%, 1.8%, 2.5%, 3.0%

Viscosity: 2,000 to > 75,000 cPs

Molecular Weight: 2 million Daltons

pH: 6.8-7.6

Osmolarity: 200 - 400 mOsM

Particulates: ≤10,000 particles>10μ

≤1,000 particles>25μ

Bacterial Endotoxins: Not more than 0.5 EU/ml

Sterility: Sterile using method described in USP<71>

**Intended use**

Viscoelastic SH is indicated for use as an ophthalmic surgical aid in the anterior segment in patients undergoing cataract extraction and intraocular lens (IOL) implantation. Viscoelastic SH is intended to create and maintain space to manipulate intraocular tissue in the anterior segment during cataract extraction and intraocular lens implantation.

Ophthalmic viscosurgical devices serve to maintain a deep anterior chamber during anterior segment surgery allowing reduced trauma to the corneal endothelium and surrounding ocular tissues. The viscoelastic properties help to push back the vitreous face and prevent formation of a flat chamber postoperatively.

**Contraindications**

At present there are no known Contraindications to the use of Viscoelastic when used as recommended; care should be used in patients with hypersensitivity to any components in this material.

**Precautions**

- Precautions normally associated with anterior segment surgical procedures should be observed.
- Postoperative increases in intraocular pressure have been reported with sodium hyaluronate products. The IOP should be carefully monitored and appropriate therapy instituted if significant increases occur.
- It is recommended that Viscoelastic SH be removed from the eye as is possible by thoroughly irrigating and aspirating with a balanced salt solution at the close of surgery.

**Warnings**

- Do not reuse cannula.
- Use only if material is clear.
- Avoid trapping air bubbles
- Use only if container is undamaged. Contents are only sterile if the pouch is undamaged.
- Do not resterilize the syringe.
- Do not use after expiration date shown on syringe as sterility cannot be assured after the expiration date.

**Adverse reactions**

- Viscoelastic SH is tolerated after injection into human eyes during ophthalmic surgical procedures. As with most viscoelastics, a transient rise in intraocular pressure has been reported in some cases.
- Postoperative inflammatory reactions such as hypopyon and iritis have been reported with the use of ophthalmic viscoelastics, as well as incidents of corneal edema and corneal decompensation. Their relationship to the use of Viscoelastic SH has not been established.

**Applications**

Cataract Surgery – IOL Implantation

Use sterile opening technique when removing contents of the peel pouches. In a sterile environment, the pouch is peeled open to attain the sterile syringe. A cannula is used to slowly and carefully inject a sufficient amount of Viscoelastic into the anterior chamber. The injection may be performed before or after removal of the crystalline lens. Additional viscoelastic will be injected as needed during anterior segment surgery to fully maintain anterior chamber depth, to manipulate intraocular tissue, or to replace any viscoelastic lost during the surgical procedure. Viscoelastic SH may also be used to coat surgical instruments and the intraocular lens prior to implantation.

**How supplied**

Viscoelastic SH is filled into a sterile, single use, pre lubricated glass syringe with Luer lock. A gray butyl stopper is inserted into the syringe, followed by a polycarbonate plunger rod. The assembled syringe is sealed into a peel pouch, steam sterilized and packaged into a box. The complete assembled box is shrink wrapped for protection. Store at approximately 25 °Celsius. Avoid exposure to direct light. Available in 1.0ml-2.5ml, accompanied with a sterile cannula.

**Directions for Use**

**FOR INTRAOCULAR USE. FOR SINGLE USE ONLY.**

1. Peel open pouch from both syringe and cannula under aseptic conditions.
2. Remove rubber tip cap from syringe tip.
3. Thread the cannula onto the syringe by holding the Luer Lock adapter between your thumb and forefinger. Twist the cannula in a continuous motion using the cartridge as a wrench with the other hand. Twist until the cannula hub has traveled the full length of the threads and is firmly seated.
4. Visually inspect that the cannula hub has traveled the full length of the Luer Lock adapter.
5. Remove the plastic cartridge from the cannula in a straight motion, being sure not to twist or unscrew the cannula while removing the cartridge.
6. Purge the remaining air from the system by holding the syringe barrel with one hand and gently depressing the plunger rod with the other until Viscoelastic SH appears at the cannula tip.

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**CAUTION**

Federal (USA) law restricts this device to sale by or on the order of a physician.

Symbol	Description
	Warning/Caution
	Do not reuse
	Batch Code
	Reference Number
	Sterilization by moist heat
	Use until
	See the instructions before use
	Don't use when packing damaged
	Manufacturer
	EC Representative
	Keep from Sunlight
	Keep dry
	Do not re-sterilize
ISO 15523-1	Symbol Reference Standard

**Manufactured By:**

EyeKon Medical, Inc.  
2451 Enterprise Road  
Clearwater, FL 33763  
Telephone (727) 799-2212  
[www.eyeconmedical.com](http://www.eyeconmedical.com)  
email:  
[info@eyeconmedical.com](mailto:info@eyeconmedical.com)  
Made in the U.S.A.

EyeKon Medical, Ltd.  
1 Friary Temple Quay  
Bristol BS1, 6EA  
United Kingdom